

Reliable | Low bleed
Robust


BP5MS GC column

Excellent performance for your MS analysis

Trajan Scientific and Medical is a vertically integrated manufacturer of gas chromatography columns, with all aspects of column manufacture controlled at one site.

Our BP5MS GC column has been designed to provide excellent performance in the areas of column bleed, activity, selectivity and robustness.

Reliable | Low bleed
Robust



GC column | BP5MS

Bleed

Our BP5MS column has been designed with low bleed performance as a criteria. Figure 1 shows the bleed performance of the BP5MS.




Figure 1. BP5MS GC column bleed performance.

Activity

Column activity is an important attribute of a high performance GC column. Limiting the active sites on the tubing (Silanol) means that reactive compounds can be chromatographed without loss in the column. The below table displays a test mix which has several probes for activity. The increasing tailing factor shows a larger activity. As can be seen the BP5MS and the market leader are almost identical in performance with the probes highlighted.

| Test mix | | BP5MS | | Market leader | |
|------------------------|--------------------|----------------|------|----------------|------|
| Compounds | Functional test | Tailing factor | k' | Tailing factor | k' |
| 1. 1-Propionic acid | Basicity | 0.62 | 0.67 | 0.62 | 0.67 |
| 2. 1-Octene | Polarity | 1.16 | 0.97 | 1.18 | 1.13 |
| 3. n-Octane | Hydrocarbon marker | 1.18 | 1.06 | 1.22 | 1.17 |
| 4. 4-Picoline | Acidity | 1.18 | 1.11 | 1.85 | 1.49 |
| 5. n-Nonane | Hydrocarbon marker | 1.92 | 1.44 | 1.21 | 1.69 |
| 6. Trimethyl phosphate | Acidity | 1.18 | 1.64 | 1.89 | 1.79 |
| 7. 1,2-Pentanediol | Silanol | 1.21 | 1.83 | 1.13 | 1.86 |
| 8. n-Propylbenzene | Hydrocarbon marker | 1.20 | 1.95 | 1.19 | 1.98 |
| 9. 1-Heptanol | Silanol | 1.21 | 2.00 | 1.20 | 2.03 |
| 10. 3-Octanone | Polarity | 1.22 | 2.11 | 1.20 | 2.13 |
| 11. n-Decane | Hydrocarbon marker | 1.19 | 2.17 | 1.18 | 2.19 |

Selectivity

Our BP5MS column has been designed to provide the same selectivity as the market leading 5MS column. Figure 2 displays the peak retention times for an EPA 8270 test mix. Both column retention time data are plotted against each other allowing a visual indication of selectivity. The mathematical correlation of 0.9996 between the two columns demonstrates the selectivity is statistically the same.




Figure 2. Selectivity of the BP5MS column compared to a market leading 5MS column.

Robustness

The robustness of the BP5MS GC column was determined by extracting 20 human serum samples (using the derivatization reagents methoxyamine hydrochloride and BSTFA+1% TMCS) and looking at the performance of specific compounds of interest to metabolomics.

- Lactic acid
- Alanine
- Valine
- Glyceric acid
- Fumaric acid
- Threonine
- Aspartic acid
- Pyroglutamic acid
- Proline
- Xylose
- Citric acid
- Glucose
- Sorbitol
- Ribose-5-phosphate
- Glucose-6-phosphate
- Cystine



Metabolite concentrations span many orders of magnitude for example glucose (mM) and glucose-6-phosphate (μM). All can be measured accurately within the $\text{CV} < 25\%$ peak area acceptable limit for biological-based analysis, as shown in figure 3 and 4.




Figure 3. Retention time reproducibility CV (%) on early eluting metabolites (e.g. lactic acid) and on highly abundant metabolites (e.g. glucose).




Figure 4. Peak area reproducibility CV (%) for 16 common metabolites – all are within the CV<25% acceptable limits for metabolomics analysis.

The same BP5MS GC column was run with 100 samples of the human serum. Figures 5, 6, 7 and 8 demonstrate the robustness of the GC column, showing the 100th sample result.




Figure 5. Example of the 100th human serum sample.




Figure 6. Separation of leucine/isoleucine in the 100th human serum sample.




Figure 7. Identification of leucine in the 100th human serum sample against the NIST 17 library.





Figure 8. Identification of isoleucine in the 100th human serum sample against the NIST 17 library.

BP5MS column

Excellent performance for your MS analysis

The BP5MS GC column has been designed to provide excellent performance across a range of applications specific for your analysis.



Visit us at www.trajanscimed.com or contact your regional Trajan representative for assistance and further information.

Trajan Scientific and Medical

Science that benefits people

Trajan is actively engaged in developing and delivering solutions that have a positive impact on human wellbeing. Our vision revolves around collaborative partnerships that improve workflows, delivering better results.

BFL-0567-G © Trajan Scientific Australia Pty Ltd 09/2019

BGB GC|LC
MS|CE

www.bgb-shop.com

- Switzerland:** BGB Analytik AG • Rohrmattstrasse 4 • 4461 Böckten • Phone +41 61 991 00 46 • Fax +41 61 991 00 25 • sales@bgb-analytik.com
BGB Analytik SA • Route de Pré-Bois 20 • 1215 Genève 15 • Phone +41 22 788 49 43 • Fax +41 22 788 49 45 • sales.fr@bgb-analytik.com
- Benelux:** BGB Analytik Benelux B.V. • Drielandendreef 42-44 • 3845 CA Harderwijk • Phone +31 341 700270 • Fax +31 341 700271 • sales.benelux@bgb-analytik.com
- France:** BGB Analytik France S.A.S. • 81 Vie de l'Etraz • 01630 St. Jean de Gonville • Phone +33 450 488567 • Fax +33 450 562378 • sales.fr@bgb-analytik.com
- Germany:** BGB Analytik Vertrieb GmbH • Mühlestraße 1 • 79539 Lörrach • Phone +49 7621 5884270 • Fax +49 7621 5884289 • sales.de@bgb-analytik.com
- Turkey:** BGB Analytik Laboratuvar Malzemeleri Ticaret A.Ş. • Cevizli Mah. 30 Ağustos Cad. No: 13 • B-Blok, Kat: 14 Daire 127 • 34846 Maltepe, İstanbul
Phone +90 2169092048 • sales.tr@bgb-analytik.com